




EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH
Steingrabenstrasse 14
DE-82211 Herrsching
Germany
Delivery Address: Am Anger 5, DE-82237 Woerthsee-Ettersschlag

Tel +49 (0) 8153 88702-00
Fax +49 (0) 8153 88702-10
www.easycap.de
info@easycap.de

Cap for Real-Behavior-EEG, 24Ch Motor Layout for Smarting



Details for Users

Ordering Information

Catalogue-Number: RBE-24-MOTOR-# (please replace # by cap size)

Cap sizes are given in cm head circumference:

- Adult caps: 54, 56, 58, 60, 62 (average male: 58, average female: 54)
- Children caps: 50 (5 years), 52 (7 years), 54 (11 years)
- Infant caps: 36 (newborn), 38, 40 (3 months), 42, 44 (7 month), 46, 48 (3 years)

Note: We supply caps with adjusted cuts for Asian head shapes. If requested, please mention it in your order.

The catalogue-number comprises the cap as described, serial number, and this document; all packed in a labelled cardboard box. For further information about accessories or consumables please visit our website or contact our local distributor.

Cap

Sublnion Cap with integrated chin belt, Black, High-Elasticity Fabric, with velcro-patch below Oz to hold the amplifier's backpack.

AsiaCut upon request.

Electrodes

All electrodes are Multitrodes with sintered Ag/AgCl sensors. They are buttoned directly into the cap (total height 3,5 mm).

All electrodes are number-labelled at the sensor end.

The cable colours correspond to the above figure.

The cables are threaded in and out of the cap and meet occipitally in a common connector.

Tornister

A Velcro-attachable Tornister for Smarting is NOT part of each cap. If you have no Tornister yet, one per Smarting needs to be ordered additionally.

Theta/Phi-Coordinates

Please find a table with Theta/Phi-Coordinates of all electrode sites at the end of this file.

Table of Coordinates

Pin-Number	Channel-Number	Name	Theta	Phi
11	1	Fp1	-90	-72
9	2	Fp2	90	72
7	3	Fz	45	90
5	4	F7	-90	-36
3	5	F8	90	36
1	6	FC1	-31	-46
2	7	FC2	31	46
4	8	Cz	0	0
6	9	C3	-45	0
8	10	C4	45	0
10	11	T7	-90	0
12	12	T8	90	0
14	13	CPz	22	-90
16	14	CP1	-31	46
18	15	CP2	31	-46
20	16	CP5	-69	21
22	17	CP6	69	-21
24	18	TP9	-113	18
26	19	TP10	113	-18
25	20	Pz	45	-90
23	21	P3	-60	51
21	22	P4	60	-51
19	23	O1	-90	72
17	24	O2	90	-72
15	CMS	CMS	23	90
13	DRL	DRL	67	90

These values are standardized to a Theta of 90° for the plane through Fpz, T7, T8, Oz.

The signs follow this convention:

